

NEUROLOGICAL:

Status: Normal / Abnormal: ? Paralyzed Vec. Cooled to 34.5
 Sedation: Morphine / Alfentanil / Fentanyl / Propofol / Midazolam
 Paralyzed: No / Yes

SEPSIS: (trend ↑ or ↓)

WCC:
 Temp: 34.5
 Platelets:
 +ve culture: Sputum / BAL / Urine / Blood Date: -
 Organism: Antibiotics: Fluclax / Nexid
 Day: 2

PLAN: (ward round)

am	pm	evening
- Continue to cool - Full inotropic support - Enoximone to 5mg/kg - Keep tvc on catheter - ↓ Fenc/Vac to 1ml/hr - Try to feed. SF		

COMMENTS:

10⁰⁰ WR with IDW
 for PD catheter to decompress abdo study (20ml/hr)

- NOT TO DRAW FREELY
 MAX DRAW 30ml/hr

For 2 x 20ml COLLOIDS - one stat
 - second over 2 hrs.

enoximone doubled

1700 20ml/hr ascates drained off.
 Colloid repl. c FFP.

Chyler
 SF